camlog

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3 Kontakt

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English

The following descriptions contain detailed instructions on cleaning, disinfection, and sterilization of the instruments and prosthetic components of the CAMLOG® Implant System. Detailed information on the choice of suitable implants, prosthetic components, treatment planning and the use of CAMLOG® implants is available in the user information and on the CAMLOG Biotechnologies AG website at www.camlog.com. Your local CAMLOG representative will be happy to provide you with information.

Important information: <p>If not otherwise specified in the instruction manual, reusable CAMLOG® products may be reprocessed as long as they are maintained in working condition according to the instruction manuals and processing procedures.</p>
<p>CAMLOG® products intended for single use must not be reused because safe preparation and/or functional safety cannot be ensured.</p>

1 Instruments

The instruments used in the CAMLOG® Implant System are supplied non-sterile unless they are explicitly marked as sterile. They must be cleaned, disinfected, and sterilized before first use and every use.

1.1 Cleaning and Disinfection

Thorough cleaning and disinfection is indispensable for effective sterilization. The following information refers to a manual preparation procedure with a combined cleaning and disinfecting agent. When choosing a combined cleaning and disinfecting agent, one should make sure that:

- it is suitable for cleaning and disinfecting dental instruments
- it is suitable for ultrasonic cleaning (no foam development)
- it has a proven efficacy for disinfection
- it is compatible with the materials of the products to be cleaned and disinfected
- it does not contain aldehyde (otherwise blood, secretions, tissue remains, etc. may stick)

Caution! The procedure below must be observed exactly.

- Initial disinfection:** Immediately after use immerse all instruments in a combined cleaning and disinfecting agent to serve as own security and to prevent the contaminants from drying. Always adhere to the manufacturer’s instructions regarding concentration and reaction time of the combined cleaning and disinfecting agent.
- Disassembling:** Completely disassemble all instruments consisting of several parts (see 1.3 Disassembling and Assembling). Remove the depth stops from the drills.
- Preliminary cleaning:** Remove coarse contaminants from the instruments within 2 hours after use. For this purpose use a disposable cloth, running water, and a soft brush (no metal bristles). Remove residues in areas with difficult access using an appropriate instrument.
- Combined cleaning and disinfection:** Place the instruments for cleaning and disinfection in the combined cleaning/disinfecting solution for the scheduled application time, making sure they are completely covered. The instruments are not to come into contact with one another. In case of remaining debris use an ultrasonic unit. Brush off the instruments with a soft brush (no metal bristles) so as to completely remove residue. The cooling channels in the internally cooled instruments and the drill extension must be rinsed out five times with the combined cleaning/disinfecting solution using the cleaning cannula (Art. No. J5002.0020) and a syringe. Remove concealed tissue remains with the cleaning needle (Art. No. J5002.0012). The driver for screw implant, long, manual/ratchet (Art. No. J5300.0009), the adapter, ISO-shaft for contra-angle/ratchet (Art. No. J5002.0010) and the driver for ball attachment manual/ratchet (Art. No. J5300.0011) have a tapped blind hole. Remove the remains contained in the tapped blind hole with an appropriate instrument. The concentrations and application times specified by the manufacturer of the combined cleaning/disinfecting agent must be observed.
- Rinsing:** Remove the instruments, rinse them thoroughly under deionized water with a low bacterial count and low endotoxin count (such as Aqua Purificata). Take particular care in to the areas with limited access such as channels or grooves.
- Drying:** Dry the instruments. Use a loose lint free disposable cloth and compressed dry air, free of oil and with a low bacterial count. We recommend the use of a sterile filter.
- Inspection:** Inspect the instruments for signs of corrosion, surface damage, chipping or contamination. Remove damaged instruments from use. Following inspection, if the instruments are visually contaminated, reclean an rerinse.
- Assembly:** Reassemble all disassembled instruments (see 1.3 Disassembly and Assembly).
- Packaging:** Pack instruments for sterilization promptly. We recommend using a disposable sterilization package. You must ensure that the sterilization package is suitable for steam sterilization (constant temperature of at least 141 °C/286 °F, sufficient vapor permeability) and that the products are adequately protected against mechanical damage.

1.2 Sterilization

Permissible steam sterilization procedures are fractionated vacuum procedures (with sufficient product drying). Other sterilization procedures are not allowed. Care should be taken that:

- the maximum steam sterilization temperature is 138 °C/280 °F
- the sterilization time (exposure time at sterilization temperature) is at least 20 minutes at 121 °C/250 °F or at least 5 minutes at 134 °C/273 °F.
- A minimum drying time of 30 minutes is recommended for each of the cycles above.

The use of a less effective gravitation procedure must be supported by the validation of the sterilizer and the validation of the control process. Care should be taken with surgery trays that they do not touch the walls of the steam sterilizer. High local temperatures (over 160 °C/320 °F) lead to deformation of plastics.

Attention: All non-sterile packaged CAMLOG® products must not be sterilized in the CAMLOG original packaging!
--

1.3 Disassembly and Assembly

The following instruments must be disassembled prior to cleaning and disinfection:

- Cardanic driver for screw implants** (Art. No. J5300.0010)

Disassemble the cardanic driver in three parts by loosening the clamping screw with a screw driver hex. After reassembling the clamping screw must be tightened by hand with a screw driver hex.

- Driver for bar abutment** (3.3 mm, 3.8 mm, 4.3 mm: Art. No. J5300.0020; 5,0 mm, 6,0 mm: Art. No. J5300.0025)

To disassemble the driver for bar abutment remove the bolt by turning it counter-clockwise. Reassemble by turning clockwise.

- Driver for impression post and healing cap for bar abutment** (3.3 mm, 3.8 mm, 4.3 mm: Art. No. J5300.0027; 5,0 mm, 6,0 mm: Art. No. J5300.0028)

To disassemble the driver for impression post and healing cap for bar abutment remove the bolt by turning it counter-clockwise. Reassemble by turning clockwise.

- Adapter for screw implants, long** (3.3 mm: Art. No. J5302.3310; 3.8 mm: Art. No. J5302.3810)

To disassemble the adapter for screw implants, remove the bolt by turning it counter-clockwise. Reassemble by turning clockwise.

- Torque wrench** (Art. No. J5320.1030)
For disassembling, maintenance and reassembling of the torque wrench follow the “Instruction Manual – CAMLOG® Torque wrench”, Art. No. J8000.0044.

2 Prosthetic Components

The prosthetic components of the CAMLOG® Implant System are supplied non-sterile unless they are explicitly marked as sterile. They must be used one time only and solely on one patient. They must be cleaned and disinfected before each use on a patient and after each use on a patient, e.g. for shipping to the dental laboratory. We recommend an additional sterilization.

Attention: Plastic components of the CAMLOG® Implant System must not be sterilized (only cleaned and disinfected), except for the temporary abutment made out of PEEK. The temporary abutment made out of PEEK can be sterilized.
--

Prosthetic components connected with a nonsterilizable material cannot be sterilized (just cleaned and disinfected).

2.1 Cleaning and Disinfection

Thorough cleaning and disinfection is indispensable for effective sterilization. The following information refers to a manual preparation procedure with a combined cleaning and disinfecting agent.

When choosing a combined cleaning and disinfecting agent, one should make sure that:

- it is suitable for cleaning and disinfecting dental-prosthetic components
- it is suitable for ultrasonic cleaning (no foam development)
- it has a proven efficacy for disinfection
- it is compatible with the materials of the products to be cleaned
- it does not contain aldehyde (otherwise blood, secretions, tissue remains, etc. may stick)

Caution! The procedure below must be observed exactly.

- Initial disinfection:** Immediately before and after use immerse all prosthetic components in a combined cleaning and disinfecting agent to serve as own security and to prevent the contaminants from drying. Always adhere to the manufacturer’s instructions regarding concentration and reaction time of the combined cleaning and disinfecting agent.
- Combined cleaning and disinfection:** Place the prosthetic components for cleaning and disinfection in the combined cleaning/disinfecting solution for the scheduled application time making sure they are completely covered. The prosthetic components are not to come into contact with one another. In case of persistent contamination use an ultrasonic unit. Brush off the prosthetic components with a soft brush (no metal bristles) so as to completely remove residues. Remove all debris contained in the areas with limited access. Observe the concentration and application times specified by the manufacturer of the combined cleaning/disinfecting agent.
- Rinsing:** Remove the prosthetic components, rinse them thoroughly under deionized water with a low bacterial count and low endotoxin count (such as Aqua Purificata). Take particular care in rinsing the area that have limited access.
- Drying:** Dry the prosthetic components. Use a loose lint free disposable cloth and compressed dry air, free of oil dry and aseptic. We recommend the use of a sterile filter.
- Inspection:** Inspect the prosthetic components for signs of corrosion, surface damage, chipping or contamination. Damaged prosthetic components cannot be used. Following inspection, if the components are visually contaminated, reclean an rerinse.
- Packaging:** Pack prosthetic components for sterilization promptly. We recommend using disposable sterilization package. You must ensure that the sterilization package is suitable for steam sterilization (constant temperature of at least 141 °C/286 °F, sufficient vapor permeability) and that the products are adequately protected against mechanical damage.

2.2 Sterilization

Permissible steam sterilization procedures are fractionated vacuum procedure (with sufficient product drying) and gravitation procedure. Other sterilization procedures are not allowed. Care should be taken that:

- the maximum steam sterilization temperature is 138 °C/280 °F.
- the sterilization time (exposure time at sterilization temperature) is at least 20 minutes at 121 °C/250 °F (gravity procedure as well as fractionated vacuum procedure) or at least 5 minutes at 134 °C/273 °F (gravity procedure as well as fractionated vacuum procedure).
- A minimum drying time of 30 minutes is recommended for each of the cycles above.

Attention: All non-sterile packaged CAMLOG® products must not be sterilized in the CAMLOG original packaging!
--

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